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Courriel : [contact@in-surg.fr](mailto:contact@in-surg.fr)

[www.in-surg.fr](http://www.in-surg.fr)

**REUSABLE SURGICAL INSTRUMENTS**
**1. INTENDED USE**

The IN'SURG instruments are reusable surgical instruments of operating room quality provided "non-sterile." They are intended without being connected to an active medical device to perform a surgical procedure (non-exhaustive list) such as: drilling, sawing, scraping, scraping, tightening, retracting, spreading, or attaching biological tissues. IN'SURG surgical instruments are intended to be reused after being subjected to appropriate cleaning, disinfection, and sterilization procedures.

**2. CONTRAINDICATIONS**

There are no known contraindications associated with the use of IN'SURG surgical instruments. However, they should only be used by a trained, qualified healthcare professional who adheres to the recommendations outlined in this manual.

**5. INSTRUCTIONS FOR TREATMENT OF REUSABLE SURGICAL INSTRUMENTS**

This section provides end users, in accordance with standard EN ISO 17664, with instructions for cleaning, disinfecting, sterilizing, and maintaining reusable IN'SURG instruments. These instructions are intended to assist surgeons and medical personnel in effectively and safely processing medical devices.

**3. WARNINGS**

- Before first use and before each subsequent use, thoroughly clean, disinfect, and sterilize all instruments and check their operation.
- Inform the manufacturer and/or the competent authority of any serious incident related to the IN'SURG surgical instrument.
- The instrument should only be used for the purpose for which it was designed and validated for the corresponding type of surgery.
- Instruments should not be used if the user or healthcare professional is informed that the patient has intolerances to the materials they contain.

**4. MATERIALS**

The main materials used in the manufacture of IN'SURG surgical instruments are stainless steel according to standard NF EN 10088-3 and titanium.

<b>General recommendations</b>	<ul style="list-style-type: none"> <li>• Prohibit the use of sodium hypochlorite (bleach) to prevent corrosion pitting on the instruments.</li> <li>• Always follow the usage and dosage instructions provided by the detergent product manufacturer.</li> <li>• Check the surface condition and functionality of the instrument before each use.</li> <li>• It is recommended to reprocess the instruments as soon as possible after use.</li> </ul>
<b>Processing limit</b>	<p>Due to the design of the product and the materials used, no limit regarding the maximum number of treatment cycles that can be performed can be specified. The lifespan of medical devices is determined by their function and careful handling.</p>
<b>Initial treatment / Pre-disinfection</b>	<ul style="list-style-type: none"> <li>• Prepare a pre-disinfection solution in accordance with the disinfectant product manufacturer's recommendations.</li> <li>• Immerse the instruments (all IN'SURG instruments are submersible), in the open position, immediately after use.</li> <li>• Maintain immersion for the duration specified by the disinfectant product manufacturer.</li> <li>• Carefully brush, with a nylon brush, the surfaces that require it, starting with those accessible in the folded state, then those accessible in the open state until all visible dirt is removed. Pay attention to hard-to-reach spaces (e.g., joints).</li> <li>• Rinse thoroughly with running water according to the disinfectant product manufacturer's recommendations.</li> </ul> <p><b>Comments:</b>  <i>Recommended product: Detergent-disinfectant without aldehyde or quaternary ammonium, bactericidal, fungicidal, yeast-killing, tuberculocidal, and effective against enveloped viruses. Follow the instructions and adhere to the concentrations and immersion times provided by the manufacturer (excessive concentration or duration can lead to corrosion or protein fixation).</i></p>
<b>Transportation</b>	<p>The transportation of the instruments from the operating room to the reprocessing area must be done in closed containers to avoid contaminating the environment.</p>
<b>Automatic Cleaning &amp; Disinfection (Washer-Disinfector)</b>	<p>Place the open instruments on the support provided in the washer-disinfector and start the cleaning cycle.</p> <ol style="list-style-type: none"> <li>1 - Pre-wash with cold water (network or softened) for a minimum of 2 minutes.</li> <li>2 - Wash at 50/55 °C (according to the technical data sheet of the detergent used), for a minimum of 5 minutes with a universal alkaline, enzymatic detergent containing surfactants in softened water.</li> <li>3 - Rinse with cold or hot softened water for a minimum of 2 minutes.</li> <li>4 - Thermal disinfection at 90°C (minimum 5 minutes) or 93 °C (minimum 3 minutes) with osmosis water.</li> <li>5 - Drying according to the manufacturer's recommendations for the washer used.</li> </ol> <p><b>Comments:</b></p> <ul style="list-style-type: none"> <li>• Follow the instructions (temperature, contact time, concentration) provided by the detergent solution manufacturer.</li> <li>• Verify the cycle's compliance with the parameters established during the annual operational qualification.</li> <li>• Follow the washer-disinfector instructions and verify the success criteria after each cycle has been completed, as indicated by the manufacturer.</li> <li>• Use only washer-disinfectors compliant with EN ISO 15883, maintained and validated regularly.</li> </ul>
<b>Functional control and testing</b>	<ul style="list-style-type: none"> <li>• Perform a visual inspection to verify the cleanliness and integrity of the instruments, their assembly, and their proper functioning.</li> <li>• If necessary, repeat the reprocessing process until the instrument is visibly clean.</li> <li>• Verify the functionality of each instrument.</li> </ul>
<b>Packaging</b>	<ul style="list-style-type: none"> <li>• Place the instruments in a container or packaging for steam sterilization in accordance with ISO 11607-1 and EN 868 standards.</li> <li>• Sterilize in an autoclave (steam):             <ul style="list-style-type: none"> <li>◦ Instruments with closed rack at the first rack notch.</li> <li>◦ Articulated instruments in the closed position.</li> <li>◦ Disassembled instruments in the assembled state.</li> </ul> </li> </ul>

<b>Sterilization</b>	<p>Steam sterilization of instruments in accordance with the national requirements in force (134° for 18 minutes). The sterilization cycle is detailed as follows:</p> <ul style="list-style-type: none"> <li>• Pretreatment step (3 successive phases of vacuum followed by steam injection to prepare the instrumentation for sterilization).</li> <li>• Vacuum before pressure and temperature rise.</li> <li>• Sterilization plateau: maintaining a minimum temperature of 134°C. Plateau holding time 18 min (French regulation)</li> <li>• Drying vacuum and maintaining this vacuum for about 20 min.</li> <li>• Injection of filtered air to return to ambient pressure before opening the door.</li> </ul> <p><b>Comments:</b></p> <ul style="list-style-type: none"> <li>• <i>The sterilization process must be validated according to standard NF EN 285 and compliant with standard NF EN ISO 17665.</i></li> </ul>
<b>Storage after reprocessing</b>	<p>Store the sterilized instruments in a dry, clean, and dust-free environment in the sterilization packaging.</p> <p><b>Comments:</b></p> <p><i>Check packaging and medical devices before use (package integrity, absence of moisture, and expiration date). In case of damage, a complete reprocessing must be carried out.</i></p>
<b>Reprocessing validation information</b>	<p>The following control instructions, materials, and machines were used for the validation of IN'SURG instruments:</p> <p>Manual disinfectant agent: Sekusept aktiv (Ecolab)  Machine cleaning agent: Neodisher MediClean forte Dr. Weigert (alkaline)  Disinfectant washer: Miele &amp; Cie. KG Type: G 7836 CD Mobile injector E440/2</p> <p><b>Comments:</b></p> <p><i>If the machines and chemicals listed above are not used, it is the responsibility of the user to validate their process according to the manufacturer's instructions for the products used. The user must ensure that the treatment method employed, including resources, materials, and personnel, is appropriate and enables compliance with applicable requirements. The state of the art and national laws require following validated procedures.</i></p>

## 6. DISPOSAL

The instruments must be disposed of in accordance with current national and local laws and regulations.

## 7. STORAGE AND TRANSPORT INSTRUCTION

There are no specific requirements for the transport and storage of reusable surgical instruments IN'SURG. However, it is recommended to store the instruments away from heat and humidity.

## 8. RETURN CONDITIONS

Before returning defective products, it is necessary to complete the entire reprocessing process (see Chapter 5).

Instruments sent for demonstration purposes may be handled, displayed, inspected, but must not be used or sterilized. They should be stored in their original packaging for return.

## 9. WARRANTY CONDITIONS

IN'SURG defines warranty periods for its instruments, with the end of the warranty year engraved on them. The warranty period varies from 2 to 10 years depending on the range. However, damages caused by the user, natural wear and tear are not covered by the warranty. Any unauthorized maintenance by IN'SURG cancels the warranty.

Any liability is excluded if the defects or their consequences result from handling or modification of the product by the customer or a third party not authorized by the company IN'SURG.

## 10. ADDITIONAL INSTRUCTIONS

**10.1.** IN'SURG surgical instruments are submersible and washable in a washer-disinfector.

- It is recommended to follow the soaking time and dosage recommended by the manufacturer of the pre-disinfection product used.
- IN'SURG surgical instruments are also compatible with manual washing provided that protocols are followed (wearing personal protective equipment, preparing the washing bath, rinsing, drying, etc.).
- The use of metal brushes, scouring pads, and abrasive products is prohibited.

**10.2.** When recomposing surgical sets:

- Ensure the cleanliness of the instruments.
- Ensure that the joints of articulated instruments are flexible. If necessary, apply a suitable lubricant to the joints that require it, manipulate the joint, and remove excess oil with a compress.

- Check the integrity and functionality of each instrument before packaging them.

**10.3.** During individual packaging:

- Check the integrity and functionality of the instruments to be packaged.
- Protect the ends of fragile instruments.

Use packaging appropriate to the size of the instrument to be packed, in accordance with current recommendations.

## 11. EXPLANATION OF SYMBOLS

	Catalog reference		Date of manufacture
	Batch code		Manufacturer
	Numéro de série		Non-sterile
	Consult the instructions for use		Fears moisture (Keep dry)
	CE marking		Protect from heat and radioactive sources
	Medical devices		Unique device identifier
	Can be sterilised in a steam steriliser (autoclave) at the specified temperature.		

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